



MARINA MEDICAL

The New Wave in OBGYN

TAKE CONTROL

Of Your Vaginal Repair Procedures


The
TAKE CONTROL SYSTEM
Includes:

- Martin's Arm Table Assistant
- 360° Rotatable Table Clamp
- VIPA™ - Fiberoptic Pelvic Assistant



- Marina Vaginal Cuff Suture Platform

*New, Elongated,
Illuminated!*



TAKE CONTROL with Martin's Arm...

...A Cost Effective, Valuable O.R. Assistant

- Attaches to the O.R. Table
- Steadier than the Human Hand
- Eliminates Shake & Fatigue
- Adjusts to Any Height or Angle
- Quick Adjust & Lock Control

110-652
or
110-655

254-620

254-516

110-652
or
110-655

Martin's Arm Table Assistant Clamp provides complete rotation whether O.R. table track is horizontal or vertical.

110-655-020
and
110-655-030

110-655
plus
110-660


Martin's Arm Table Assistant Set

Description

Catalog #

Martin's Arm with Clamp - 22 in. x 10 in., L-Stand (extended), 45 in. height	110-652
Martin's Arm Table Assistant, 33 in. height (Table Clamp #110-660 not included)	110-655
Martin's Arm Retractor Table Clamp	110-660
Mini Mayo Tray for Martin's Arm	110-655-020
Tray Holder for Martin's Arm	110-655-030
Steiner Speculum 1.1 Kg/2.5 lb; 14cm x 3cm /5.5 in. x 1.25 in.	252-055
Breisky Retractor - 140mm x 35mm/5.5 in. x 1.3 in.	254-516
Breisky-Navratil Retractor - 170mm x 25mm/6.75 in. x 1 in.	254-620

110-660




TAKE CONTROL with VIPA™


Veronikis Intra-Fiberoptic Pelvic Assistant

The Most Exciting New Instrument Available for:

- Sacrocolpopexy
- Paravaginal Repair
- Urethropexy
- Rectovaginal Septum Endometriosis



*Now Featuring New Elongated,
Illuminated
Suture Platforms*



New Suture Platforms are 12 centimeters in length allowing illuminated examination of vaginal tissue and a continuous suturing platform for the entire length of the vagina.

Illumination Helps:

- Find the Apical Scar
- Pinpoint suturing on steady platform
- Avoid suturing thinned out tissue

Marina Vaginal Cuff Suture Platforms:

The Set
Actual Size
Diameters



19 mm



25 mm



31 mm

VIPA™ Suture Platforms

Description

- 0.75 in / 19 mm Diameter, 12 cm length
- 1.00 in / 25 mm Diameter, 12cm length
- 1.25 in / 31 mm Diameter, 12cm length

Catalog #

- 350-403
- 350-404
- 350-405

TAKE CONTROL

with other

WIPA™ Tips



Description

Catalog #

VIPA Urethropexy Replacement Tip	350-374
VIPA Vaginal Cuff Suture Platform (small)	350-376
VIPA Vaginal Cuff Suture Platform (large)	350-378
VIPA Rectal Probe	350-380

Marina Vaginal Stents

For hand held use in the O.R. during vaginal surgery. Marina Vaginal Stents are 10 inches long and come in a set of three different diameters:

- 18 mm (3/4 in.)
- 25 mm (1 in.)
- 31 mm (1 1/4 in.)

Custom sized diameters are available upon request.



Marina Lucite Vaginal Stents

Description

Catalog #

0.75 in / 19 mm Diameter, 10" length	350-382
1.00 in / 25 mm Diameter, 10" length	350-384
1.25 in / 31 mm Diameter, 10" length	350-386

MARINA MEDICAL

955 Shotgun Road, Sunrise, Florida 33326 USA

(954) 924-4418 • Fax (954) 924-4419 • (800) 697-1119 • Fax (800) 748-2089 • Se habla español • www.marinamedical.com